

Technical data

Dimensions	60(B) x 200(H) x 370(D) mm
Weight	4,0kg
Segment lenght	65/70/75/80/85/90/95mm

Power supply

Voltage	100-240 VAC
Frequency	47-63Hz
Power consumption (sealing)	200VA max
Power consumption (on standby)	30VA max
Power fuse	3,15A/slow

Sealing head &

RF generator

Nominal output power	50W for load res. 50 Ohm
Sealing frequency	40,68 MHz±5kHz
Standard sealing time	1-3 seconds
Sealing time limitation	7,5 sec ± 1 sec
Sealing jaws gap (opened)	5 ±0,2mm
Sealing jaws gap (closed)	0,07mm

Environmental

Temperature Working conditions	10°C - +40°C
Temperature Storage conditions	-20°C - +50°C